

TA Prime

Food Texture Analyzer



MADE IN THAILAND

The Stealth TA Prime Series

Table Top Food Texture Analyzer are designed to test a wide spectrum of materials including low force application for a wide range of food types, including baked goods, cereals, confectionaries, snacks, dairy, fruit, vegetables, gelatins, meat, poultry, fish, pasta and even pet food & Packaging etc.

Ordering Information



Tech-Quality Co.,Ltd.



Call : +662-3286878



2/13 Sukapiban 2 Rd.,
Pravet,Pravet, Bangkok 10250
Thailand



+662-3286818



www.tech-quality.com




info@tech-quality.com

@techquality



In food texture testing, standard tests such as compression, tension and flexure are used to measure Brittleness, Chewability, Chewiness, Elasticity, Fracturability, Gumminess, Hardness1, Hardness2, Rubberiness, Springiness, Stringiness, Adhesive, Adhesiveness, Consistency, Crispness, Crumb Extensibility, Crunchiness, Extensibility, Gel Strength, Hardness work done, Index of viscosity, Ripeness, Rupture Strength, Spreadability, Toughness, Work of Cutting, Work of Failure. Work of Penetration, Work of Shear



■ Simple to use Test Setup

Limit Testing Break Testing
Loop Testing Preload Setup
Elongation Auto Return,
Over Load Protection,
Auto Save after test

■ Multi-control mode

- Constant speed
- Constant Displacement
- Multi-step control

■ Applications

- Food & Packaging, etc.
- Bulk Analysis
- Tension / Compression
- Extrusion / Back Extrusion
- Puncture & Penetration
- Shear
- Bend & Break
- Peel & Seal

■ Specification

TA Prime Machine Capacity : 1KN
DC Motor Drive set
Crosshead Speed Range : 15-800mm/min
Displacement Resolution : 0.00001mm
Crosshead Travel : 0.01 - 300 mm
Calibration Tension & Compression 10 point
Load Cell Accuracy : +/- 0.1% FS
/ ISO17025 – 1 Class 0.5
Speed Accuracy : 2%
Data Transfer : Rate to PC 1 kHz.
Operating Temperature : 5° to 35°C (41° to 95°F)
Communication port : USB
Supply Voltage : 220Vac ±10% 50-60 Hz 3Amp
High Accuracy : Ball Screw Type
Weight Approx : 30 kg
Warranty : 1 Year




Blackhole software

Multi Graph Displays



Excel , PDF User design reports



Display	Variable	Unit	Decimal
Display 1	Force	N	2
Display 2	Hardness1	gf	3
Display 3	Crispness	-	3
Display 4	Stringiness	s	6
Display 5	Adhesiveness	N/s	3

Multi Display & Decimal Point

Shape: Circular

Sample Name: Crocodile

Unit: mm

Parameter(mm): Width: 0.0, Thickness: 0.0, ID: 0.0, OD: 0.0

Gauge Length: 50.000, Diameter: 10.000, Span length: 0.0

Peel Set: StartPoint: 4.00, EndPoint: 25.00

Temperature: 25

Area: 78.540 mm

Auto Selete Tab Number: 100

Sample Information

Step	Direct	Control Break Type	Speed:mm/min	Target:Unit
1	Tensile	PreLoad	1000:mm/min	1:N
2	Tensile	Load	50:mm/min	15:N
3	Tensile	Stroke	300:mm/min	29:mm
4	Tensile	Stress	20:mm/min	10:kgf/mm.2
5	Compress	Load	99:mm/min	1:kgf

Multi-stage testing with Bravo software

Header	Value
Address	2/13 Sukapiban 2 Rd.,Pravet , Bangkok 10250 Thailand.
Tel	+662 3286878 (Auto)
Website	www.tech-quality.com
ID LINE	@techquality

Statistic Reports

Auto Zero before Test: ☒

Auto Save after Test Complete: ☒


Auto Return to home after Test: ☒

Auto Return to Limit after Test: ☐


Setting: Return Speed: 900.0 mm/min, Load Protect: 80 % 16.00 kgf

Sensor: ☐ Stroke ☐ Rub Extensometer ☐ Extensometer


Auto test setup & load protect



“Formula design by user”




“Selectable result”



“Results can be color-coded for immediate pass/fail recognition”



Library



Standard Library



Probe Guide



YouTube Guide


Food Texture Analyzer



Wafer



Noodle



Rice



Bread



**Tomato
Ketchup**



**Ice Cream
Cone**



Crispy Rice



**Cabbage and
Chinese Cabbage**



Hot Dog



Crispy Cookies



Mantou



Cereals